COLLECTION AND SHIPMENT OF FIBROBLAST SAMPLES FOR IFAR REGISTRY THE ROCKEFELLER UNIVERSITY HOSPITAL

Laboratory of Genome Maintenance Agata Smogorzewska, MD, PhD, Phone: 212-327-7850

PURPOSE:

This sample is being obtained **for research purposes.** Any genetic testing that might be done on this sample has to be confirmed in a clinical laboratory. As the sample is being obtained for research purposes, we cannot say how long it will take to discover findings. Sometimes it can take months or years before we are able to provide any information to the referring physician.

CONSENT:

No sample should be collected or shipped unless a valid consent is obtained. For participants already enrolled in IFAR, consent may already be on file. For questions about consent status, please contact our study coordinator (contact information below).

NOTIFICATION:

Please notify our study coordinator (contact information below) at least one week in advance of the scheduled skin biopsy procedure OR 2 days in advance for the shipment of cultured fibroblasts.

	SKIN BIOPSY	CULTURED FIBROBLASTS
COLLECTION	New fibroblast primary cultures are set up	Provide two T25 flasks of cultured fibroblasts
	at Integrated Genetics, which are then	at 80% confluency (at time of shipment), at
	returned to our lab. Collection and shipping	room temperature.
	materials will be provided by IFAR. Collect	
	one 3 to 4 mm punch and place directly into	If fibroblasts have already been frozen please
	the sterile transport medium provided.	send one frozen vial containing at least one
		million viable cells on dry ice.
TIMING	Samples should be collected/shipped MONDAY-THURSDAY. Samples should NOT be	
	collected/shipped on Friday, Saturday, or Sunday.	
PACKING	It is extremely important to carefully pack the specimens, providing ample protection against	
	breakage. Use an insulated secondary container and include absorbent material AND sturdy	
	outer packaging as an overpack. FedEx air bills will be provided by IFAR, but please	
	complete Sender's Info before shipping.	
TEMPERATURE	Shipped at room temperature	Cultured Fibroblasts at room temp
		Frozen Fibroblasts on dry ice
REQUISITION	The IFAR will provide you with the	Please complete the requisition form provided
	completed requisition form for Integrated	by the IFAR and include it with the cells sent
	Genetics. This must accompany the biopsy	directly to the IFAR.
	being sent	
SHIPPING	FED EX PRIORITY OVERNIGHT	FED EX PRIORITY OVERNIGHT
	Integrated Genetics/Genzyme	The Rockefeller University
	2000 Vivigen Way	Laboratory of Genome Maintenance
	Santa Fe, NM 87505	Attn: Frank Lach
	800-848-4436	1230 York Avenue Box 182
	Integrated Genetics sends cultured cells to	New York, NY 10065
	our laboratory for all analyses	212-327-8862

QUESTIONS/CONTACT INFORMATION:

For questions about the research study, contact our genetic counselor: Jennifer Kennedy at 212-327-8612. For questions about specimen retrieval/shipment, please contact: Frank Lach at 212-327-8862.