


ATP/NADH Coupled ATPase Assay Protocol with Synergy Neo Multi-Mode Reader



Monitor NADH Conc. by either of these:

1. Absorbance at 340nm measured by *Spectrometer* (Extinction coefficient of NADH: $6220 \text{ M}^{-1}\text{cm}^{-1}$)
2. Fluorescent emission at 445nm (340 excitation) measured by NEO Plate Reader

Stock solution to make:

1M HEPES, pH 8.0
50mM NADH (Roche, 4°C. Freshly made or store stock at -80°C)
0.2M ATP (Sigma, -20°C. pH to 7.0 with KOH and store at -80°C)

5mL Reaction Buffer:

250uL 1M HEPES, pH 8.0
30uL Pyruvate Kinase (PK), 10mg/mL stock (Roche, 4°C)
32uL Lactate Dehydrogenase (LDH), 5mg/mL stock (Roche, 4°C)
9.3mg phosphoenolpyruvate (PEP), powder, (Roche, 4°C)
15uL (30uL) 50mM NADH
(50uL 0.1M MgCl₂, added according to your needs)
4.658 ddH₂O

Final Conc.

50mM
60mg/mL
32mg/mL
9mM
0.15mM (0.30 mM)
(1mM)

OPTIONAL: DTT or buffers at different pH

To Measure ATPase activities with NEO Plate Reader:

- 1) Generate a standard curve of NADH fluorescence vs. NADH concentration at time 0 with only NADH in buffer. (Make sure NADH conc. used in the reaction buffer falls in the linear range.)
- 2) To measure ATPase activities:
 - a. Sample: mix 0.5mL reaction with desired amount of your ATPase
(For PCAT1, 0.066mg/mL protein was used in the reaction to measure basal activities). Control/background: mix 0.5mL reaction with H₂O instead of ATPase.
- 3) Aliquot 30uL reaction mix or control to each well. (Make sure to measure enough duplicates to eliminate errors due to pipetting or air bubbles.)
- 4) Add 1-2 uL of desired amount of ATP to each well by 384 well multi-channel pipette.
- 5) Run ATPase coupled protocol saved in the machine with desired temperature, duration, and frequency of measurement. (Get training for the plate reader from another lab member).

- 6) Export the measurements and plot both control and real measurement in excel.
- 7) Calculate background NADH oxidation rate with the control by calculated the slope of fluorescent decrease over time.
- 8) Calculate ATPase rate from the slope of fluorescent decrease over time and substrate away the background NADH oxidation rate.
- 9) Convert the ATPase rate to nmol/mg/min:
 - a. Calculate the rate NADH disappearance (mM/min) from the standard curved measured at 1).
 - b. Convert mM/min to umol/mL/min
 - c. NADH disappearance= ATPase rate, so divide the rate by amount of ATPase used (mg/mL) to get the rate with the unit of nmol/mg/min

- 10) If ATPase in the proteoliposomes were used in the measurements, divide the rate by the ratio of ATPase facing outside if it's known. Otherwise multiple the rate by 2 fold assuming half of the ATPases are inside, and half are outside of the liposome.