

## PCR Reaction

## PCR Conditions

### **For Phusion DNA Polymerase (Thermo Fisher)**

	<b>Standard</b>	
5X Phusion HF buffer *	- 10 ul	98°C – 1 min
10mM dNTP	- 1 ul	98°C – 15 s
Primers/each (10µM)	- 2.5 ul	55°C – 20 s
cDNA <b>OR</b>	- 0.5 ul	72°C –*
Genomic DNA <b>OR</b>	- 0.5 ul	72°C – 4 min
Plasmid DNA	- 1.0 ul	4°C - for ever
(1:50 dilution of miniprep DNA ~5ng/ul)		
Water	- 32.5-33 ul	
Phusion Turbo S	- 0.5 ul	

50ul

\*7.5mM MgSO<sub>4</sub> in 5X Phusion HF buffer

\*extension time: 30s per 1kb

### **For Deep Vent DNA Polymerase (NEB)**

	<b>Standard</b>	
10X ThermoPol Buffer	- 5 ul	95°C – 2 min
10mM dNTP	- 1 ul	95°C – 20 s
Primers/each (10µM)	- 2.5 ul	55°C – 20 s
cDNA <b>OR</b>	- 0.5 ul	72°C –*
Genomic DNA <b>OR</b>	- 0.5 ul	72°C – 6 min
Plasmid DNA	- 1.0 ul	4°C - for ever
(1:50 dilution of miniprep DNA ~5ng/ul)		
MgSO <sub>4</sub> (25mM)	- 3ul	
Water	- 34.5-35 ul	
Deep Vent	- 1 ul	

50ul

\*extension time: 1min per 1kb

After PCR, run the products on 1% Agarose gel. If product is correct size, then cut off the band and use NucleoSpin Gel Clean-up Kit to purify DNA. For higher DNA yields,

\* Use hot ddH<sub>2</sub>O (~60°C) and incubate for ~4 min before elution.

## Restriction Enzyme Digests( for E.coli cloning)

Digest with appropriate NEB Enzymes and Buffers (see Tools&Resources on the NEB website and then select Double Digest Finder to learn the information of specific restriction enzyme ).

- \* Using larger reaction volume (~35-40 ul) is helpful
- \* For digest time, 2-3 hr is enough in most cases.

Gel-purify cut PCR product using NucleoSpin Gel Clean-up Kit.

Elute with 15-30ul ddH<sub>2</sub>O.

Measure concentration and label tube.

### **Vectors**

- To digest a vector for cloning, first digest with selected NEB enzymes, then CIP (removes phosphates from vector for better cloning) treat.
- Example: pNEH cut with EcoRI and Sall

DNA	1-2ug plasmid DNA (~5ul of a miniprep) or 20.25ul cleaned up PCR rxn
Enzyme 1 (EcoRI)	1ul
Enzyme 2 (Sall)	1ul
10X Buffer	2.5ul
100X Bsa	0.25 ul
Water	to 40ul
<hr/>	
	40ul

- Digest for 2hrs at 37°C, then for vectors only add 0.2ul CIP (NEB) for 10min at 37°C.
- Immediately gel-purify the cut vector using a NucleoSpin Gel Clean-up Kit
- Elute with 15-30ul ddH<sub>2</sub>O.
- Measure concentration and label tube.

## Ligation and transformation ( for E.coli cloning)

### Ligation

- T4 ligase and T4 ligation buffer (10X)
- Mix the following contents in order and then react for ~12hr at 16°C:

○ Vector	*
○ Gene Insert	*
○ Water	to 20ul
○ 10X Ligation Buffer	2ul
○ T4 Ligase	1ul
-----	
	20ul

\*Note that vector and gene insert should have a 1:4-1:10 molar ratio (i.e. 4 pieces of insert for each piece of vector). Always use 10-20ng of vector.

### Example:

- 4Kb Vector and 1Kb Insert have a ratio of 1:4.
- If you use 15ng of vector [ $15 \times 0.25 = 3.75$ ]
- Then you would need to use 15ng of insert [ $3.75 \times 4 = 15$ ]

### Transformation

- Thaw E.coli competent cells (XL2-Blue or TOP10) on ice.
- Pipet 5-20 ul ligation reaction directly into cells and mix by tapping gently. DO NOT mix by pipetting up and down.
- Incubate on ice for 30min.
- Heat shock cells in 42°C water bath for 45-60s.
- Return tube to ice.
- Add 350ul LB Media to cells and shake at 37°C for 1hr.
- Spread 50-200 ul from each transformation vial on separate, labeled LB plate with appropriate antibiotic and incubate at 37°C overnight.
- Select colonies and analyze by PCR or sequencing.

Sam Zhao (Modified from Steve Long's protocol)  
April 2016[Type here]

\*Note that XL2-Blue cells are chloramphenicol resistant. If you are using chloramphenicol resistant vectors, such as pACYC Duet, please select TOP10 or other competent cells.

## DNA Sequencing

- Visit website <https://clims3.genewiz.com/default.aspx>
- Create a username and password
- To submit a sample click on “Creat Sequencing Order”
- Fill in appropriate information for each sequencing item
- Click Submit
- Label the top of all microtubes with corresponding numbers provided
- Sample amounts:
  - Plasmid DNA                      5-10 ul (50-70ng/ ul)
  - Primer                                5-10 ul (10 uM)
- After receiving the sequencing results, carefully check both sequencing and tracing file. Normally, the available sequencing length is 600-700 nt. If needed, using software such as SeqMan or Vector NTI to assembly the forward and reverse DNA fragments.