

DH10Bac Cells Protocol

*Use autoclaves glass/plastics

*Keep everything under a flame

Buffers:

TFB1 (200mL)

30mM KOAC

100mM RbCl

10mM CaCl₂

50mM MnCl₂

15% glycerol

pH to 5.8 by 1M HAC (acetic acid)

filter

store at 4C

TFB2 (50mL)

10mM MOPS

75mM CaCl₂

10mM RbCl

15% glycerol

pH to 6.5 1/ 2N HCl

filter

store at 4C

Plates → LB + Bactoagar (to make 100mL aka 5 plates: 2g LB pellets + 0.75g Bactoagar)
100ul tetracycline + 50ul kanamycin after agar cools, pour plates

“Transformation” → Scrape some cells off frozen pellet (from ‘Cells for Cloning’ rack-DH10Bac box-Blue top tube)

Add cells to 1mL of SOC in 1.5mL tube, plunging pipette tip up and down

Rotate for 1 hour at 37C

Plate 100ul and 200ul of that on 2 plates

Incubate o/n at 37C

Next day:

Put plate in cold room first thing in the morning (wrap in plastic wrap)

Set up overnight: 5mL LB broth + 5ul tetracycline + 2.5ul kanamycin + 3 colonies in a 50mL conical tube. Can do up to three tubes to ensure growth in at least one.

Shake overnight at 37C

Next day:

Prepare 40mL LB broth in a 500mL glass flask (2x)

Add 40ul tetracycline + 20ul kanamycin + 200ul from overnight growth

Grow to OD 0.6-0.8 (should take 3-4 hrs, sometimes 5hrs)

Chill cells on ice for 10min, transfer to 50mL Falcon tubes

Spin down @ 3,500rpm @ 4C for 10min

Resuspend cells gently with 16mL TFB1 Buffer in each of the tubes. DO NOT vortex. Leave on ice for 5min

Spin down 3,500rpm @ 4C for 10min

Resuspend cells with 1.6mL TFB2 Buffer (in each of the 2 tubes)

Aliquot cells by 50ul in 1.5mL tubes

Freeze in liquid nitrogen

Store in -80C

Final step → Ask someone to check the cells and compare growth to last time. DO NOT distribute cells to the lab until they're checked.

DH10Bac Cells Protocol (basic version)

6/3/99, from Blobel Lab, Ning Zheng based on protocol from Zhiping

Basically, 40 mL LB is used for one batch of preparation of one type of cells. 2X40 ml could be used too. 40 ml LB prep gives 1.6mL competent cells, which will be aliquoted into 50 uL. Each aliquot will be used for one transformation reaction.

Two buffers will be used: TFB1, TFB2.

Other stuff needed: liquid **N₂**, ice,

- Sterile TFB1(RbC1 transformation), total 200 ml

30 mM KOAC

100 mM RbC1

10 mM CaC1₂

50 mM MnC1₂

15 % glycerol(v/v)

adjust pH to 5.8 by 1M acetic acid (~80 uL) sterilize by filtration, store at 4°C indefinitely

- Sterile TFB2, total 50 mL 10mM MOPS

75 mM CaCl₂ 10 mM RbCl

0.105 g

0.551 g (CaC1₂.2H₂O)

0.060 g

7.5 ml(100%)

15% glycerol(v/v)

adjust pH to 6.5 with 2N HCl (~10 uL)

0.588 g

2.42 g

0.296 g (CaC1₂.2H₂O) 1.98 g (MnC1₂. 4H₂O)

30 ml(100%)

sterilize by filtration, store at 4 °C indefinitely

Inoculate (3 colonies) into 40 ml culture, grow at 37°C until OD₆₀₀ = 0.6 – 0.8

Chill the cells on ice for 10 min, transfer to 50 mL falcon tubes (prechill and mark microcentrifuge tubes)

Spin down at 3500 rpm at 4°C for 10 min.

Resuspend the cells gently by 16 mL (2/5 v/v) TFB1 buffer and incubate 5 min on ice.

Spin down at 3500 rpm at 4°C for 10 min. Resuspend the cells by 1.6 mL (1/25) TFB2 buffer.

Aliquot the cells into 50 μ L in microcentrifuge tubes. Freeze the cells in liquid N₂. Store at -80 °C.