

Purification of His-TEV (S 219V)-Arg

Cells: BL21(DE3) RIL containing pRK793 (glycerol stock in -80°C freezer)

Buffers to prepare:

Lysis buffer: 50mM Tris, pH 7.0, 150mM NaCl, 10% glycerol

Wash buffer: lysis buffer + 10mM imidazole

Pre-elution buffer: lysis buffer + 50mM imidazole (adjust to pH 7.0)

Elution buffer: lysis buffer + 200mM imidazole (adjust to pH 7.0)

Protein expression conditions

1. Grow cells at 37°C in LB containing 100ug/ml ampicillin and 30ug/ml chloramphenicol.
2. Induce cells when OD600 reaches 0.5-0.6 with 1mM IPTG. Lower the temperature to 30°C and grow cells for another 4 hours.
3. Harvest cells with PBS, TBS or lysis buffer at 4000g, 15min, 4°C and store cells in -80°C.

Protein prep and purification

1. Re-suspend cells in lysis buffer: 50mM Tris, pH 7.0, 500mM NaCl, 10% glycerol (20ml/L culture).
2. Lyse cells in a metal beaker with a sonicator on ice water. Set the amplitude at 50% with 0.5 second ON time/ 1.5 second OFF time and 5 minute total sonication time.
3. Centrifuge the lysed cells at 70,000g for 40 min.
4. Apply supernatant to pre-equilibrated Talon Co²⁺ resin column.
5. Wash extensively with >20CV of wash buffer. Check the flow-through with protein assay dye and wash until no color change.
6. Elute with pre-elution buffer (50mM imidazole) until no color change of protein assay dye.
7. Elute proteins with elution buffer.
8. Dialyze the elution against 20mM Tris, pH 7.0, 150NaCl, 10% glycerol and 5mM DTT O/N. (<5mM imidazole)
9. Spin down precipitated proteins, determine protein concentration, aliquot and freeze at -80°C.

Optional: Purification with AKTA

1. Concentrate TEV to 10-30mL, remove salt and load several times on cation ion exchange column Source 15Q. TEV is eluted between 150mM -350mM NaCl. (Buffer A: 20mM Tris, 7.0, 50mM NaCl, 1mM DTT; Buffer B: 20mM Tris, 7.0, 1M NaCl, 1mM DTT)