

## Protocol for high efficiency E. coli chemically competent cells

### 1. SOB medium (250ml)

Tryptone	5g	2%
Yeast extract	1.25g	0.5%
NaCl	0.146g	10mM
KCl	0.0466g	2.5mM
MgCl <sub>2</sub> •6H <sub>2</sub> O	0.508g	10mM
MgSO <sub>4</sub>	0.301g	10mM
<hr/>		
Total volume		250ml

### 2. TB buffer (200ml)

10mM	PIPES	0.604g
15mM	CaCl <sub>2</sub> •2H <sub>2</sub> O	0.441g
250mM	KCl	3.7275g
	• Add H <sub>2</sub> O to 150ml	
	• Adjust pH to 6.7 by 5M KOH	
55mM	MnCl <sub>2</sub> •4H <sub>2</sub> O	2.18g
H <sub>2</sub> O	to 200ml	

## Method

1. Pick a single colony from a LB/Tet (or other requisite antibiotic) plate and inoculate to 3ml SOB medium + 12.5µg/ml tetracycline (or other antibiotic) and culture at 37°C for 6-7 hours
2. Transfer the 3ml SOB/culture to 250ml SOB and shake at 18°C for 20-24 hours
3. When OD<sub>600</sub> reach 0.4-0.6\*, put the culture vessel to ice-water bath for 10 minutes
4. Chill down centrifuge for 15min, then harvest cells by 2500g centrifugation for 15min at 4°C
5. Pour off the medium and resuspend pellet with 40mL cold TB buffer on ice
6. Harvest cells by 2500g centrifugation for 15min at 4°C
7. Pour off the medium and resuspend cells with 10mL TB buffer on ice and add 0.75mL DMSO
8. dispense aliquot of the cells into pre-chilled 1.5ml tube and freeze the cells with nitrogen liquid, store them in -80°C